GB

Instructions for Use
POLIDIOX PLUS
Polydioxanone (PDO) Surgical Suture
SYNTHETIC, ABSORBABLE, MONOFILAMENT, USP/EP

DESCRIPTION
POLIDIOX PLUS (polydioxanone) monofilament synthetic absorbable suture is prepared from the polyester, poly(p-dioxanone). The empirical molecular formula of the polymer is (C4H6O3). POLIDIOX PLUS sutures are dyed violet (solvent violet 13 D&C No.2 C.I. 60725) to enhance visibility in tissue. Polydioxanone polymer has been found to be nonantigenic, nonpyrogenic and elicits only a slight tissue reaction during absorption. POLIDIOX PLUS suture provides an antibacterial effect by means of a broad spectrum antibacterial agent Chlorhexidine Diacetate (CHA) content not more than 60µg/m.

INDICATIONS
POLIDIOX PLUS monofilament synthetic absorbable sutures are indicated for use in all types of soft tissue approximation, including use in paediatric cardiovascular tissue where growth is expected to occur and ophthalmic surgery. POLIDIOX PLUS suture is not indicated in adult cardiovascular tissue, microsurgery and neural tissue. These sutures are particularly useful where the combination of an absorbable suture and extended wound support up the four weeks is desirable.

ACTIONS
POLIDIOX PLUS synthetic absorbable suture has been formulated to provide wound support through an extended healing period. The results of implantation studies of POLIDIOX PLUS monofilament suture in animals indicate that approximately 70% of its original strength remains two weeks after implantation. POLIDIOX PLUS suture provides tissue supporting during 60 days post implantation. Absorption is essentially complete within six months (180-240 days). POLIDIOX PLUS suture has been shown to inhibit microbial adherence and colonization of the suture by E.coli, Klebsiella pneumoniae, Staphylococcus epidermidis, Staphylococcus aureus and its Methicillin resistant strains used in-vitro and in-vivo studies.

CONTRAINDICATIONS
These sutures, being absorbable, are not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and are not to be used in conjunction with prosthetic devices, i.e., heart valves or synthetic grafts. POLIDIOX PLUS suture should not be used in patients with known allergic reactions to Chlorhexidine Diacetate (CHA).

WARNINGS
The safety and effectiveness of POLIDIOX PLUS (polydioxanone) sutures have not been established in neural tissue, adult cardiovascular tissue, or for use in microsurgery. Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

PRECAUTIONS
Adequate knot security requires the accepted surgical technique of flat and square ties with additional throws as warranted by surgical circumstances and the experience of the surgeon. As with any suture, care should be taken to avoid damage to needle and suture. Avoid the crushing or crimping application of surgical instruments, such as needle holders and forceps. Conjunctival and vaginal mucosal sutures remaining in place for extended periods may be associated with localized irritation. Subcuticular sutures should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption. Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

ADVERSE REACTIONS
Due to prolonged suture absorption, some irritation and bleeding has been observed in the conjunctiva and mild irritation has been observed in the vaginal mucosa.

HOW SUPPLIED
POLIDIOX PLUS sutures are available as sterile, monofilament, dyed (violet) strands in sizes 5/0 to 2 (metric sizes 1-5), and sterile, in a variety of lengths, as non-needed or attached to stainless steel needles of varying types and sizes. POLIDIOX PLUS sutures are available in one, two and three dozen boxes. POLIDIOX PLUS suture is supplied sterile. POLIDIOX PLUS suture is for single use only.

STORAGE
Store below 25 °C and keep away from sunlight. Protect from humidity. Do not use after expiry date.

SYMBOLS USED ON LABELLING

Do not reuse

Dyed, Absorbable, Monofilament

Do not resterilize

Catalogue Number

Do not use if package is damaged

Store below 25°C

Manufacturer

Keep away from sunlight

Date of Manufacture, Year

Protect from humidity

Expiry date, Year-Month

Recyclable pack

Sterile EO: Ethylene oxide

Attention, See instruction for use

Batch Number

IFU-PXP-rev-04-19-01-2023 Issue date: 11.09.2012

"EASSI(The European Association of Surgical Suture Industry) has developed a system of symbols which is designed to describe various suture product characteristics in an intuitive, pictorial manner. The use of symbols is permitted by the Medical Device Directive (MDD 93/42/EEC) and enables companies to provide information to the customer without having to provide multilingual translations."

BOZ TIBBİ MALZEME SANAYİ VE TİCARET A.Ş.
Sağlık Mahallesi Sağlık 1 Sk. No:335 06410 Çankaya/ANKARA/TÜRKİYE
Tel : +90 (312) 254 03 40 (5 hat) Faks : +90 (312) 254 03 50
web:www.boztibbi.com e-mail:boz@boztibbi.com

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